


**Improved Resolution of 1,4-Dioxane-d8 from the Solvent Peak when the Starting Oven Temp. is Reduced to 60 °C (Top) from 70 °C (Bottom)**



**Peaks**

1. 1,4-Dioxane-d8 (IS)
2. N-Nitrosodimethylamine
3. Pyridine

**Column** Rxi-5Sil MS, 20 m, 0.15 mm ID, 0.15 µm (cat.# 43816)  
**Sample** 8270 MegaMix (cat.# 31850)  
 8270 Benzidines mix (cat.# 31852)  
 Benzoic acid (cat.# 31879)  
 Revised B/N surrogate mix (cat.# 31888)  
 Acid surrogate mix (4/89 SOW) (cat.# 31063)  
 Revised SV internal standard mix (cat.# 31886)

**Diluent:** Methylene chloride  
**Conc.:** 20 µg/mL (IS/SS 20 µg/mL)  
**Injection**  
 Inj. Vol.: 1 µL split (split ratio 20:1)  
 Liner: Topaz 4 mm single taper w/wool (cat.# 23303)  
 Inj. Temp.: 275 °C

**Split Vent**  
 Flow Rate: 12 mL/min  
**Carrier Gas** He, constant flow  
 Flow Rate: 0.72 mL/min

**Detector** MS  
**Mode:** Scan  
**Transfer Line**  
 Temp.: 280 °C  
**Analyzer Type:** Quadrupole  
**Source Temp.:** 330 °C  
**Quad Temp.:** 180 °C  
**Electron Energy:** 70 eV

**Solvent Delay**  
 Time: 1.3 min  
**Tune Type:** DFTPP  
**Ionization Mode:** EI  
**Scan Range:** 39-550 amu  
**Scan Rate:** 9.8 scans/sec

**Instrument** Agilent 7890B GC & 5977A MSD  
**Notes** Analyzed using a 120 V oven equipped with the GC Accelerator kit (cat.# 23849).

**Oven Programs**

**Top:** 60 °C (hold 0.7 min) to 285 °C at 39.8 °C/min to 305 °C at 4.3 °C/min to 330 °C at 28.5 °C/min (hold 3.5 min)  
**Bottom:** 70 °C (hold 0.7 min) to 285 °C at 39.8 °C/min to 305 °C at 4.3 °C/min to 320 °C at 28.5 °C/min (hold 0.7 min)