


# Bile Acids with Isomer Separation in Human Plasma on Raptor Inert ARC-18 by LC-MS/MS



Peaks	tr (min)	Precursor Ion	Product Ion
1. GUDCA	1.32	448.40	74.15
2. TUDCA	1.40	498.10	80.05
3. GCA	1.78	464.10	74.15
4. TCA	1.80	514.20	80.05
5. UDCA	2.63	391.50	391.50
6. GCDCA	2.66	448.40	74.15
7. TCDC	2.66	498.10	80.05
8. GDCA	2.88	448.40	74.15
9. TDCA	2.88	498.10	80.05
10. CA	2.92	407.20	407.20
11. HDCA	2.99	391.50	391.50
12. TLCA	3.87	482.10	80.05
13. GLCA	4.00	432.20	74.15
14. CDCA	4.51	391.50	391.50
15. DCA	4.68	391.50	391.50
16. LCA	6.59	375.40	375.40
17. DHLCA	6.68	373.20	373.20

**Column:** Raptor Inert ARC-18 (cat.# 9314A12-T)  
 Dimensions: 100 mm x 2.1 mm ID  
 Particle Size: 2.7 µm  
 Pore Size: 90 Å  
 Guard Column: Raptor Inert ARC-18 EXP guard column cartridge 5 mm, 2.1 mm ID, 2.7 µm (cat.# 9314A0252-T)  
 Temp.: 50 °C  
**Standard/Sample**  
 Diluent: 60:40 Water:mobil phase B  
 Conc.: 5 µM  
 Inj. Vol.: 5 µL  
**Mobile Phase**  
 A: 5 mM Ammonium acetate in water, pH unadjusted  
 B: Methanol:acetonitrile (v/v, 50:50)

Time (min)	Flow (mL/min)	%A	%B
0.00	0.5	60	40
6.00	0.5	30	70
7.00	0.5	20	80
7.50	0.8*	0	100
8.10	0.8*	0	100
8.20	0.5	60	40
9.50	0.5	60	40

Max Pressure: 365 bar

**Detector:** Shimadzu LCMS-8045 in ESI- mode  
**Instrument:** Shimadzu Nexera X2  
**Sample Preparation:** For control samples, a 90 µL aliquot of 2x charcoal-stripped plasma (K<sub>2</sub>EDTA) was added to a microcentrifuge tube and spiked with 10 µL of calibrator/QC material and vortexed. Ten microliters of internal standards was added and vortexed. Samples were protein precipitated using 400 µL of ice-cold acetonitrile. After vortexing and centrifugation at 4200 rpm for 15 minutes, the supernatant was transferred to a glass test tube and dried down under nitrogen. All samples were reconstituted in 200 µL of 60:40 water:mobil phase B. The sample was transferred to a clean 2 mL screw-thread vial (cat.# 21143) with a glass insert (cat.# 21776) and capped with short-cap, screw-vial closures (cat.# 24498).  
**Notes:** \*The flow rate was increased to 0.8 mL/min to more thoroughly flush phospholipids from the analytical column, thereby reducing matrix effects.

The flow was diverted to waste before 1 minute and after 7 minutes to protect the mass spectrometer.