


# Analysis of 13 Alkaloids Found in Psychedelic Mushrooms Using Force Biphenyl by HPLC-UV

- Rapid 9-minute cycle time
- Full resolution of all 13 analytes



Peaks	tr (min)
1. N,N-desmethyl Psilocybin (Norbaeocystin)	1.710
2. N-desmethyl Psilocybin (Baeocystin)	2.258
3. 5-hydroxy-N-Methyltryptamine (5-hydroxy NMT)	2.636
4. Aeruginascin	2.721
5. Psilocybin	2.851
6. 5-hydroxy-N,N-Dimethyltryptamine (Bufotenine)	3.031
7. 4-hydroxy-N-Methyltryptamine (Norpsilocin)	3.231
8. 4-hydroxy-N,N-Dimethyltryptamine (Psilocin)	3.698
9. 4-hydroxy-N,N,N-Trimethyltryptamine (4-hydroxy TMT)	3.942
10. 5-methoxy-N-Methyltryptamine (5-methoxy NMT)	4.424
11. 5-methoxy-a-Methyltryptamine (5-methoxy AMT)	4.675
12. 4-acetoxy-N,N-Dimethyltryptamine (Psilacetin)	4.796
13. N,N-Dimethyltryptamine N-oxide (N,N-DMT N-oxide)	4.916

**Column** Force Biphenyl (cat.# 962931E)  
**Dimensions:** 100 mm x 3.0 mm ID  
**Particle Size:** 3 µm  
**Pore Size:** 100 Å  
**Guard Column:** Force Biphenyl EXP Guard Cartridge 5 mm, 3.0 mm ID (cat.# 962950253)  
**Temp.:** 60 °C

**Standard/Sample**

**Diluent:** Water  
**Conc.:** 50 ppm  
**Inj. Vol.:** 1 µL

**Mobile Phase**

**A:** Water, 10 mM ammonium formate, 0.1% formic acid  
**B:** Methanol, 0.1% formic acid

Time (min)	Flow (mL/min)	%A	%B
0.00	0.8	100	0
5.00	0.8	55	45
5.01	0.8	5	95
6.00	0.8	5	95
6.01	0.8	100	0
9.00	0.8	100	0

**Detector** UV/Vis @ 222 nm

**Flow Cell Size:** 500 nL

**Instrument** Waters ACQUITY UPLC H-Class

**Sample Preparation** Standards were prepared in a 2.0 mL amber, short-cap vial (cat.# 21142) with a 250 µL glass insert (cat.# 24518) and capped with a 9 mm short cap (cat.# 24497).

**Notes**

System and column were passivated using LC Passivation Solution (cat.# 32475). The outlined method conditions were used to perform five injections of passivation solution (5 µL) on column, diverting to waste. Five warm-up injections were completed prior to analysis.

**Storage Conditions:** Flush column using acetonitrile with 0.1% formic acid for 5 min.