


Rxi-5ms columns have the lowest bleed among all major brands of columns.



Column	Rxi-5ms, 30 m, 0.25 mm ID, 0.25 μ m (cat.# 13423)
Sample	1 ng tridecane on-column
Conc.:	1 ng tridecane on-column
Injection	split
Oven	
Oven Temp.:	350 °C
Carrier Gas	H ₂
Detector	FID
Notes	Comparison of 30 m x 0.25 mm ID, 0.25 μ m columns at 330 °C through 350 °C.