


Semivolatiles on Rxi-SVOCms

- Low phase ratio (β) is ideal for large analyte lists.
- Thick film provides higher capacity.
- High inertness ensures symmetrical phenol and benzidine peaks.



GC_EV1505

Column Standard/Sample

Rxi-SVOCms, 30 m, 0.25 mm ID, 0.50 μ m (cat.# 16638)
8270/Appendix IX kit (cat.# 31815)
Revised B/N surrogate mix (cat.# 31888)
Acid surrogate mix (cat.# 31063)
Dichloromethane
20 μ g/mL (2 ng on-column)

Diluent: Conc.: Injection

Inj. Vol.: 1 μ L split (split ratio 10:1)
Liner: Topaz 4.0 mm ID single taper inlet liner w/ wool (cat.# 23303)
Inj. Temp.: 250 $^{\circ}$ C

Oven

Oven Temp.: 40 $^{\circ}$ C (hold 1 min) to 280 $^{\circ}$ C at 12.4 $^{\circ}$ C/min to 330 $^{\circ}$ C at 3.3 $^{\circ}$ C/min (hold 3.65 min)
Carrier Gas: He, constant flow
Flow Rate: 1.2 mL/min
Detector: MS
Mode: Scan

Scan Program:

Transfer Line Temp.:

280 $^{\circ}$ C
Analyzer Type: Quadrupole
Source Type: Inert
Drawout Plate: 6 mm ID
Source Temp.: 330 $^{\circ}$ C
Quad Temp.: 180 $^{\circ}$ C
Electron Energy: 70 eV
Solvent Delay Time: 2 min
Tune Type: DFTPP
Ionization Mode: EI
Instrument
Sample Preparation

Agilent 7890A GC & 5975C MSD

Samples were aliquoted into amber 2 mL, 9 mm short-cap, screw-thread vials (cat.# 21143) containing glass Big Mouth inserts (cat.# 21782) and sealed with 2.0 mL, 9 mm short-cap, screw-vial closures (cat.# 23842).

Group	Start Time (min)	Scan Range (amu)	Scan Rate (scans/sec)
1	2.0	50-500	3.6



© 2022 Restek Corporation. All rights reserved.
www.restek.com

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

Semivolatiles on Rxi-SVOCms

Peaks	tr (min)
1. 1,4-Dioxane-d8	3.44
2. 1,4-Dioxane	3.48
3. N-Nitrosodimethylamine	3.87
4. Pyridine	3.94
5. Ethyl methacrylate	4.50
6. 2-Picoline	4.96
7. N-Nitrosomethylethylamine	5.08
8. Methyl methanesulfonate	5.52
9. (SS) 2-Fluorophenol	5.77
10. N-Nitrosodiethylamine	6.14
11. Ethyl methanesulfonate	6.59
12. Benzaldehyde	7.19
13. (SS) Phenol-d6	7.23
14. Phenol	7.26
15. Aniline	7.36
16. Bis(2-chloroethyl) ether	7.45
17. Pentachloroethane	7.46
18. 2-Chlorophenol	7.57
19. 1,3-Dichlorobenzene	7.85
20. (IS) 1,4-Dichlorobenzene-D4	7.95
21. 1,4-Dichlorobenzene	7.98
22. Benzyl alcohol	8.15
23. 1,2-Dichlorobenzene	8.25
24. 2-Methylphenol	8.33
25. Bis(2-Chloroisopropyl)ether	8.41
26. N-Nitrosopyrrolidine	8.59
27. 4-Methylphenol	8.61
28. 3-Methylphenol	8.61
29. N-Nitrosodi-N-propylamine	8.64
30. Acetophenone	8.65
31. 4-Nitrosomorpholine	8.66
32. o-Toluidine	8.71
33. Hexachloroethane	8.87
34. (SS) Nitrobenzene-D5	8.93
35. Nitrobenzene	8.96
36. N-Nitrosopiperidine	9.23
37. Isophorone	9.40
38. 2-Nitrophenol	9.56
39. 2,4-Dimethylphenol	9.61
40. Benzoic acid	9.69
41. Bis(2-chloroethoxy)methane	9.80
42. 2,4-Dichlorophenol	10.00
43. 1,2,4-Trichlorobenzene	10.18
44. Phentermine	10.18
45. (IS) Naphthalene-D8	10.30
46. Naphthalene	10.34
47. 4-Chloroaniline	10.42
48. 2,6-Dichlorophenol	10.44
49. Hexachloropropene	10.51
50. Hexachlorobutadiene	10.58
51. Caprolactam	11.05
52. N-Nitroso-N-butylamine	11.08
53. 1,4-Phenylenediamine	11.09
54. 4-Chloro-3-methylphenol	11.34
55. Isosafrole I	11.50
56. 2-Methylnaphthalene	11.69
57. 1-Methylnaphthalene	11.88
58. Hexachlorocyclopentadiene	12.00
59. 1,2,4,5-Tetrachlorobenzene	12.01
60. Isosafrole II	12.08
61. 2,4,6-Trichlorophenol	12.23
62. 2,4,5-Trichlorophenol	12.29
63. (SS) 2-Fluorobiphenyl	12.41
64. Safrole	12.53
65. Biphenyl	12.61
66. 2-Chloronaphthalene	12.65
67. 1-Chloronaphthalene	12.69
68. Diphenyl ether	12.82
69. 2-Nitroaniline	12.83
70. 1,4-Naphthoquinone	12.98
71. 1,2-Dinitrobenzene	13.06
72. 1,3-Dinitrobenzene	13.10

Peaks	tr (min)
73. Dimethyl phthalate	13.20
74. 2,6-Dinitrotoluene	13.32
75. 1,4-Dinitrobenzene	13.41
76. Acenaphthylene	13.47
77. 3-Nitroaniline	13.64
78. (IS) Acenaphthene-d10	13.75
79. Acenaphthene	13.82
80. 2,4-Dinitrophenol	13.85
81. 4-Nitrophenol	13.95
82. Pentachlorobenzene	14.08
83. 2,4-Dinitrotoluene	14.12
84. Dibenzofuran	14.16
85. 1-Naphthalamine	14.31
86. 2,3,5,6-Tetrachlorophenol	14.31
87. 2,3,4,6-Tetrachlorophenol	14.39
88. 2-Naphthalamine	14.46
89. Diethyl Phthalate	14.61
90. Fluorene	14.85
91. 4-Chlorophenyl phenyl ether	14.85
92. 2-Methyl-5-nitroaniline	14.85
93. 4-Nitroaniline	14.85
94. 4,6-Dinitro-2-methylphenol	14.92
95. N-Nitrosodiphenylamine	15.07
96. N,N-Diphenylhydrazine	15.16
97. (SS) 2,4,6-Tribromophenol	15.32
98. 1,3,5-Trinitrobenzene	15.57
99. Phenacetin	15.67
100. Diallate	15.72
101. 4-Bromophenyl phenyl ether	15.82
102. Hexachlorobenzene	15.94
103. Atrazine	16.14
104. Pentachlorophenol	16.32
105. 4-Aminobiphenyl	16.33
106. Pentachloronitrobenzene	16.36
107. Propyzamide	16.46
108. (IS) Phenanthrene-D10	16.72
109. Phenanthrene	16.77
110. Anthracene	16.86
111. Carbazole	17.18
112. di-n-Butyl phthalate	17.89
113. 4-Nitroquinoline 1-oxide	18.29
114. Isodrin	18.85
115. Fluoranthene	19.14
116. Benzidine	19.40
117. (SS) Pyrene-D10	19.56
118. Pyrene	19.59
119. Aramite I	19.87
120. (SS) p-Terphenyl-d14	19.91
121. Aramite II	20.02
122. Dimethylaminoazobenzene	20.19
123. 4,4'-Dichlorobenzilate	20.28
124. 3,3'-Dimethylbenzidine	20.87
125. Butyl benzyl phthalate	20.91
126. Kepone	20.99
127. Bis(2-ethylhexyl) adipate	21.09
128. 2-(Acetylamino)fluorene	21.41
129. 3,3'-Dichlorobenzidine	22.11
130. Benz[a]anthracene	22.18
131. (IS) Chrysene-D12	22.21
132. Chrysene	22.27
133. Bis(2-ethylhexyl) phthalate	22.28
134. Di-n-octyl phthalate	24.14
135. Benzo[b]fluoranthene	25.24
136. 7,12-Dimethylbenzo[a]anthracene	25.24
137. Benzo[k]fluoranthene	25.33
138. Benzo[a]pyrene	26.30
139. (IS) Perylene-D12	26.48
140. 3-Methylcholanthrene	27.50
141. Dibenz[a,j]acridine	29.58
142. Indeno[1,2,3-cd]pyrene	30.18
143. Dibenz[a,h]anthracene	30.28
144. Benzo[ghi]perylene	31.11