

Rum CarboBlack™ B

5% Carbowax® 20M 80/120 CarboBlack™ B (cat.# 80105)
2m, 1/8" OD x 2mm ID SilcoSmooth™ tubing
0.5µL on-column injection, on-column conc.: neat

Oven temp.: 65°C (hold 5 min.) to 150°C @
4°C/min.
Inj./det. temp.: 200°C / 250°C
Det.: FID
Carrier gas: nitrogen
Column flow rate: 20mL/min.
FID sensitivity: 1.28×10^{10} AFS

1. acetaldehyde
2. methanol
3. acetone
4. ethyl formate
5. ethanol
6. ethyl acetate
7. *n*-propanol
8. *sec*-butanol
9. isobutanol
10. active amyl alcohol
11. isoamyl alcohol
12. *n*-amyl alcohol

