

Butter Triglycerides Rtx[®]-65TG

30m, 0.25mm ID, 0.10 μ m Rtx[®]-65TG (cat.# 17008)
0.2 μ L cold on-column injection of 1% butterfat in isoctane.

Oven temp.: 80°C (hold 1 min.) to 240°C @ 30°C/min.
to 360°C @ 4°C/min. (hold 5 min.)
Det. temp.: 380°C
Carrier gas: hydrogen
Linear velocity: 70cm/sec.
FID sensitivity: 16 x 10⁻¹¹ AFS
SPI injector: high performance capillary insert 60°C,
300°C/min. to 400°C (hold 5 min.)

