

Antiarrhythmic Digoxin on Ultra IBD

Peak List:

1. unknown
2. unknown
3. digoxin

Sample:

Inj.: 10 μ L
Conc.: 1000 μ g/mL
Solvent: water:methanol (1:1, v/v)

Column:

Ultra IBD
Cat. #: 9175565
Dimensions: 150 x 4.6mm
Particle size: 5 μ m
Pore size: 100Å

Conditions:

Mobile phase: water with
0.1% (v/v) acetic acid:
acetonitrile (65:35, v/v)
Flow rate: 1.0mL/min.
Temp.: 27°C
Det.: UV @ 220nm

