

Sugars on Ultra Amino

Sample:

Inj.: 20 μ L
Solvent: mobile phase

Column:

Cat. #: 9107365
Dimensions: 150 x 4.6mm
Particle size: 3 μ m
Pore size: 100Å

Conditions:

Mobile phase: water:acetonitrile (25:75, v/v)
Flow: 0.8mL/min.
Temp.: 35°C
Det.: refractive index @ 35°C

Peak List:

	Conc.(mg/mL)
1. fructose	2.0
2. glucose	2.1
3. sucrose	4.0
4. maltose	4.5
5. lactose	4.4

