

Organic Acids Standard on Allure™ Organic Acids

Peak List:	Conc. (mg/mL)
1. tartaric acid	1
2. quinic acid	1
3. malic acid	1
4. citric acid	1
5. fumaric acid	0.005

Sample:

Inj.: 10 μ L standard solution
Solvent: water

Column:

Allure™ Organic Acids
Cat. #: 9165585
Dimensions: 300 x 4.6mm
Particle size: 5 μ m
Pore size: 60Å

Conditions:

Mobile phase: 100mM phosphate buffer, pH 2.5
Flow: 0.5mL/min.
Temp.: ambient
Det.: UV @ 226nm

